



## Vitrex Sterilance Lite3

### Sterile Safety lancet

Vitrex Sterilance Lite3 lancet is intended for blood glucose and other in vitro diagnostic tests. Health workers worldwide face accidental needle-stick injury while caring for patients, putting themselves at risk to bloodborne illnesses including the AIDS virus and Hepatitis.

#### **BUTTON ACTIVATION**

**USER SAFETY:** Needle is fully shielded before and after use.

**USER CONVENIENCE:** Pre-loaded, pressure activated design.

**PATIENT COMFORT:** High quality, ultra sharp tri-bevel needle or bi-bevel blade, high speed penetration.

**EFFICACY:** Variety of point configurations and penetration depths to suit most capillary blood requirements.

**PATIENT SAFETY:** Single-use only. Sterilised by gamma irradiation.



1



2



3



4



## Recommended use

1. Twist off the protective cap.
2. Select fingertip insertion site slightly off centre to minimise any discomfort and maximise blood flow. Hold the lancet firmly against site and press the release button to activate the device.
3. Gently massage the finger in the direction of the puncture site to obtain the required volume of blood.
4. Dispose the safety lancet in a suitable safety container.

## Technical Specifications

Colour Code	Specification	Penetration Depth	Article No.	Indicative use
Green	18G blade	1.8 mm	373213	For heel puncture in average-weight newborns and older children
Orange	21G needle	1.8 mm	377213	Haemoglobin (blood banks), cholesterol, blood typing, coagulation, blood gas and other diagnostic tests
Pink	21G needle	2.4 mm	378213	As the Orange lancet, but for thicker skin
Dark Blue	23G needle	1.8 mm	379213	Blood glucose testing and tests requiring larger samples
Yellow	26G needle	1.8 mm	375213	Diabetes clinics/hospital wards for blood glucose monitoring
Blue	26G needle	2.4 mm	376213	As the Yellow lancet, but for thicker skin
Rose	28G needle	1.5 mm	371213	Diabetes clinics/hospital wards for blood glucose monitoring
Purple	28G needle	1.8 mm	374213	As the Rose lancet, but for thicker skin
Beige	30G needle	1.8 mm	372213	Diabetes clinics/hospital wards for blood glucose monitoring. For minimal pain



Transit packaging      Pallet

**Packaging:**  
 Primary packaging: 100 pcs.  
 Transit packaging: 1,000 pcs.  
 Pallet: 90,000 pcs.

Shelf life 5 years