



Vitrex® Octaheel Baby and Vitrex® Steriheel Baby Sterile Cutting Safety lancets

Vitrex baby lancets are intended for high volume blood flow with minimum depth of penetration. The safety lancets are primarily applicable for newborns and preemies. However, they can also be recommended for other applications, where short penetration depth and/or minimal compression of tissue is important.

Vitrex Octaheel Baby and Vitrex Steriheel Baby lancets have the same indications and features, but Vitrex Octaheel is the newest edition, which offers improved ergonomics.



Vitrex Steriheel Baby

Vitrex Octaheel Baby

ACTIVATION BUTTON: Secured by safety tab.

NO CHOCKING HAZARD: Non-detachable safety tab.

USER SAFETY: Blade is fully shielded before and after use.

EASIER HANDLING: Anti-slip texture and one hand operation.

USER CONVENIENCE: Pre-loaded, button activated design.

PATIENT SAFETY: Single-use only. Sterilised by gamma irradiation.

PATIENT COMFORT: Ultra sharp blade with high speed cutting securing minimal trauma and fast healing.

NARROWED LANCING SLOT: Easier positioning.

1



2



3



4



5



LANCETS

Recommended use

1. Release the trigger lock mechanism by pushing down the safety tab.
2. Reference position for infant heel blood sampling.
3. Raise the baby's foot to a suitable position, place the blade slot flush against the heel site and activate the lancet.
4. Collect the blood sample with a capillary tube.
5. Dispose the safety lancet in a suitable safety container.

WARNING:

Do not use, if the safety tab has previously been activated or if the shield blister packing has been damaged.

Technical Specifications

Colour Code	Specification	Octaheel Article No.	Steriheel Article No.	Application
Green ●	0.65 mm Depth / 1.40 mm Incision Length	335313	335213	Micro-Premie
Orange ●	0.85 mm Depth / 1.75 mm Incision Length	334313	334213	Premie
Yellow ●	1.00 mm Depth / 2.50 mm Incision Length	333313	333213	Newborn
Pink ●	1.14 mm Depth / 2.80 mm Incision Length	332313	332213	Newborn Plus
Blue ●	2.00 mm Depth / 3.00 mm Incision Length	331313	331213	Toddler



Transit packaging

Pallet

Packaging:

Blister Packed
 Primary packaging: 50 pcs.
 Transit packaging: 200 pcs.
 Pallet: 23,200 pcs.

Shelf life 5 years