



## Vitrex No Additive Tubes

Vitrex No Additive Tubes do not contain any additives. The tubes are used to collect samples of body fluids (blood, urine, cerebrospinal fluid, exudates, etc.). The tubes can be used for specimen analysis, storage and transport. The range of application is: analysis of body fluids in clinical biochemistry, serology and immunology.

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Vacuseance+ Article No.	Vacuseance Article No.
Z/No Additive	2	13 x 75	No Additive	Paper	101102+ *	101102 *
Z/No Additive	4	13 x 75	No Additive	Paper	102204+ *	102204 *
Z/No Additive	6	13 x 100	No Additive	Paper	103106+ **	103106 **
Z/No Additive	9	16 x 100	No Additive	Paper	104109+ ***	104109 **

**Shelf life 18 months for Vacuseance**  
**Shelf life 12 months for Vacuseance+**

### Vacuseance+



Transit packaging

Pallet

**Packaging:**  
 Primary packaging: 50/100 pcs.  
 Transit packaging: 800/1,000 pcs.  
 \* Pallet: 54,000 pcs.  
 \*\* Pallet: 48,000 pcs.  
 \*\*\* Pallet: 38,400 pcs.

### Vacuseance



Transit packaging

Pallet

**Packaging:**  
 Primary packaging: 100 pcs.  
 Transit packaging: 1,200 pcs.  
 \* Pallet: 48,000 pcs.  
 \*\* Pallet: 38,400 pcs.