



Vitrex Sedimentation Systems

Vitrex Vacusence Sedimentation Tubes contain 3.2% sodium citrate for ESR test. The tubes are available with an open or closed sedimentation system.

Description	Draw Vol. (ml)	Size (mm)	Additive	Label	Vacusence+ Article No.	Vacusence Article No.
ESR Tube	2	13 x 75	4NC/Sodium Citrate 3.2%	Paper	111602+	111602
ESR Rack	2				-	116601 *

Shelf life 24 months for Vacusence
Shelf life 14 months for Vacusence+

Vacusence+



Transit packaging Pallet

Packaging:
Primary packaging: 100 pcs.
Transit packaging: 1,200 pcs.
Pallet: 54,000 pcs.

Vacusence



Transit packaging Pallet

Packaging:
Primary packaging: 100 pcs.
Transit packaging: 1,200 pcs.
* Pallet: 42,000 pcs.
** Pallet: 48,000 pcs.