



## Vitrex® Self-Sealing Haematocrit Tubes

Vitrex Self-Sealing Haematocrit Tubes offers a combination of safety and efficiency, as it contains triple Mylar wrap and a vented, self-sealing plug. The Mylar wrap provides enhanced safety with no compromise in accuracy. In the event of breakage, the Mylar wrapping will help contain the liquid and provide an added layer of protection against cuts from broken glass. The tubes are efficient as the user saves the time and trouble of having to seal the tubes with wax since it is pre-installed. To dispense the tube, the Vitrex Self-Sealing Tube Dispenser is needed. The tubes are available plain and heparinised. They are packed in vials of 100 pcs.

## Vitrex® Self-Sealing Tube Dispenser

Vitrex Self-Sealing Tube Dispenser is designed to dispense Vitrex Self-Sealing Haematocrit Tubes. In the need of centrifuging the tubes to separate the plasma from the red blood cells, then dispensing Vitrex Self-Sealing Haematocrit Tubes with other devices can be quite difficult. Vitrex Self-Sealing Tube Dispenser use positive displacement and therefore the process is a lot more convenient for the user. The dispenser is packed in boxes including 1 piece.

**For more information please contact us.**

