



## Chirana® 2-part disposable syringes and Chirana® 2-part disposable syringes with needles

Chirana 2-part disposable syringes without and with needle are made of polypropylene/polyethylene and rubber-free. They are pyrogen-free, latex-free, DEHP-free, PVC-free and non-toxic. The syringes are available with Luer tip. The syringe with needle is available in other needle combinations.

**USER CONVENIENCE:** Ergonomic design, smooth flow, highly transparent barrel, tight, black extended graduation, easy to read, double stop lines and guaranteed 0 ml position.

**USER SAFETY:** Individually blister packed and the needle is side-packed in a single blister package.

**COMFORT:** Designed tip for best connection, double stop lines and scratch-free for glove manipulation.

**PATIENT SAFETY:** Single-use only. Sterilised by ETO.

### Chirana 2-part disposable syringes

Vol. (ml)	Pcs./Box	Pcs./Carton	Article No.
2	100	2800	CH002L
5	100	1900	CH005L
10	100	1200	CH010L
20	80	960	CH020L

**Shelf life 5 years**

### Chirana 2-part disposable syringes with needles

Specification	Pcs./Box	Pcs./Carton	Article No.
2 ml with needle 0.60 x 30 mm (23G x 1 1/4")	100	2800	CH002L0630
5 ml with needle 0.70 x 35 mm (22G x 1 2/5")	100	1900	CH005L0735
10 ml with needle 0.80 x 40 mm (21G x 1 1/2")	100	1200	CH010L0840
20 ml with needle 0.80 x 40 mm (21G x 1 1/2")	80	960	CH020L0840

**Shelf life 5 years**



According to EN, ISO standards. Quality management certified according to EN ISO 13485:2016. CE marked.