

CHIRASORB braided

CHIRASORB braided is a mid-term absorbable, coated, synthetic, multifilament suture made of 90% glycolide and 10% L-lactide. The substances contained in the coating and suture are nonantigenic and noncollagenous.

Sterile Chirasorb braided is available in violet or natural. The product elicits a minimal acute inflammatory reaction in tissues.

Progressive loss of tensile strength and absorption of Chirasorb braided occurs by hydrolysis, when the polymer is degraded to glycolic and lactic acids, which are subsequently absorbed and metabolised by the body.



Indication: General soft tissue approximation and/or ligation in all specialties, including use in Ophthalmic surgery.

Contraindication: Chirasorb braided should not be used in tissues subject to expansion, stretching or distension, or which might require long-term mechanical support.

Sizes: 10/0 to 6 USP

Structure: Multifilament

Absorption method: Hydrolysis

Strength Retention

Profile: MID-TERM HEALING (35 days)

1 Day 100%

14 Days 75%

21 Days 50%

Mass Absorption

profile: 56 - 70 days



CHIRASORB braided violet ligatures



UK04
66

Pre-cut threads

EP	USP	Length	Cat. No
1	5/0	3 x 45 cm	LV 0012
1.5	4/0	3 x 45 cm	LV 0013
2	3/0	3 x 45 cm	LV 0014
3	2/0	3 x 45 cm	LV 0015
3.5	0	3 x 45 cm	LV 0016
4	1	3 x 45 cm	LV 0017
5	2	3 x 45 cm	LV 0018
1	5/0	6 x 45 cm	LV 0042
1.5	4/0	6 x 45 cm	LV 0043
2	3/0	6 x 45 cm	LV 0044
3	2/0	6 x 45 cm	LV 0045
3.5	0	6 x 45 cm	LV 0046
4	1	6 x 45 cm	LV 0047
5	2	6 x 45 cm	LV 0048
2	3/0	10 x 45 cm	LV 0064
3	2/0	10 x 45 cm	LV 0065
3.5	0	10 x 45 cm	LV 0066
4	1	10 x 45 cm	LV 0067
5	2	10 x 45 cm	LV 0068
3	2/0	12 x 45 cm	LV 4565
5	2	12 x 45 cm	LV 4568
1	5/0	2 x 75 cm	LV 0022
1.5	4/0	2 x 75 cm	LV 0023
2	3/0	2 x 75 cm	LV 0024
3	2/0	2 x 75 cm	LV 0025
3.5	0	2 x 75 cm	LV 0026
4	1	2 x 75 cm	LV 0027
5	2	2 x 75 cm	LV 0028
1.5	4/0	5 x 75 cm	LV 0053
2	3/0	5 x 75 cm	LV 0054
3	2/0	5 x 75 cm	LV 0055







Pre-cut threads

EP	USP	Length	Cat. No
3.5	0	5 x 75 cm	LV 0056
4	1	5 x 75 cm	LV 0057
5	2	5 x 75 cm	LV 0058
2	3/0	6 x 75 cm	LV 0059
3	2/0	6 x 75 cm	LV 0051
3.5	0	6 x 75 cm	LV 0060
4	1	6 x 75 cm	LV 4557
1.5	4/0	1 x 150 cm	LV 0083
2	3/0	1 x 150 cm	LV 0084
3	2/0	1 x 150 cm	LV 0085
3.5	0	1 x 150 cm	LV 0089
4	1	1 x 150 cm	LV 0087
5	2	1 x 150 cm	LV 0090

Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
2	3/0	1 x 250 cm	LV 1004
3	2/0	1 x 250 cm	LV 1005
3.5	0	1 x 250 cm	LV 1006
4	1	1 x 250 cm	LV 1007
5	2	1 x 250 cm	LV 1008
1.5	4/0	1 x 300 cm	LV 0003
2	3/0	1 x 300 cm	LV 0004
3	2/0	1 x 300 cm	LV 0005
3.5	0	1 x 300 cm	LV 0006
4	1	1 x 300 cm	LV 0007
5	2	1 x 300 cm	LV 0008
6	3 + 4	1 x 300 cm	LV 0009

CHIRASORB braided







Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 10 	0.7	6/0	1 x 45 cm	Violet	LV 0246	
	1	5/0	1 x 75 cm	Violet	LV 0248	
HR 13 	0.7	6/0	1 x 45 cm	Violet	LV 0251	
	1	5/0	1 x 45 cm	Violet	LV 0252	
	1.5	4/0	1 x 45 cm	Violet	LV 0253	
	2	3/0	1 x 45 cm	Violet	LV 0250	
	0.7	6/0	1 x 75 cm	Violet	LV 1251	
	1	5/0	1 x 75 cm	Violet	LV 0254	
	1.5	4/0	1 x 75 cm	Violet	LV 0255	
HR 16 	1	5/0	1 x 75 cm	Violet	LV 3265	
HR 18 	1.5	4/0	1 x 45 cm	Violet	LV 0258	
	2	3/0	1 x 45 cm	Violet	LV 0259	
	0.7	6/0	1 x 75 cm	Violet	LV 0241	
	1	5/0	1 x 75 cm	Violet	LV 0242	
	1.5	4/0	1 x 75 cm	Violet	LV 0256	
	2	3/0	1 x 75 cm	Violet	LV 0257	
	3	2/0	1 x 75 cm	Violet	LV 0260	
HR 20 	1.5	4/0	1 x 75 cm	Violet	LV 1261	
	2	3/0	1 x 75 cm	Violet	LV 1259	
HR 22 	1.5	4/0	6 x 45 cm	Violet	LV 0526	
	2	3/0	6 x 45 cm	Violet	LV 0527	
	1.5	4/0	1 x 75 cm	Violet	LV 0261	
	2	3/0	1 x 75 cm	Violet	LV 0262	
	3	2/0	1 x 75 cm	Violet	LV 0263	
	2	3/0	4 x 75 cm	Violet	LV 0322	
	3	2/0	4 x 75 cm	Violet	LV 0323	
	2	3/0	1 x 120 cm	Violet	LV 0363	







CHIRASORB braided

UK04
68

Absorbable - Chirasorb braided

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 22 A 	2	3/0	6 x 45 cm	Violet	LV 0327	Take-off needle
HR 25 	5	2	1 x 75 cm	Violet	LV 0297	
HR 26 	3	2/0	1 x 45 cm	Violet	LV 0136	
	2	3/0	1 x 75 cm	Violet	LV 0535	
	3	2/0	1 x 75 cm	Violet	LV 0536	
	3.5	0	1 x 75 cm	Violet	LV 0540	
HR 26 A 	2	3/0	6 x 45 cm	Violet	LV 0347	Take-off needle
	1.5	4/0	1 x 45 cm	Violet	LV 0245	
	2 3	3/0 2/0	6 x 45 cm 6 x 45 cm	Violet Violet	LV 0537 LV 0538	
HR 27 	1.5	4/0	1 x 75 cm	Violet	LV 0265	
	2	3/0	1 x 75 cm	Violet	LV 0266	
	3	2/0	1 x 75 cm	Violet	LV 0267	
	3.5	0	1 x 75 cm	Violet	LV 0268	
	4	1	1 x 75 cm	Violet	LV 1269	
	2 3	3/0 2/0	1 x 90 cm 1 x 90 cm	Violet Violet	LV 0366 LV 0367	
HR 27 A 	3	2/0	6 x 45 cm	Violet	LV 0338	Take-off needle
	2	3/0	4 x 75 cm	Violet	LV 0639	Take-off needle






CHIRASORB braided

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
2 x HR 27 	2	3/0	1 x 90 cm	Violet	LV 2466	
	3	2/0	1 x 90 cm	Violet	LV 2467	
	2	3/0	1 x 45 cm	Violet	LV 0169	
HR 30 	2	3/0	1 x 75 cm	Violet	LV 0269	
	3	2/0	1 x 75 cm	Violet	LV 0273	
	3.5	0	1 x 75 cm	Violet	LV 0274	
	4	1	1 x 75 cm	Violet	LV 1275	
	2	3/0	1 x 90 cm	Violet	LV 0369	
	3	2/0	1 x 90 cm	Violet	LV 0364	
	3.5	0	1 x 90 cm	Violet	LV 0371	
	4	1	1 x 90 cm	Violet	LV 0372	
HR 35 	2	3/0	1 x 75 cm	Violet	LV 0270	
	3	2/0	1 x 75 cm	Violet	LV 0271	
	3.5	0	1 x 75 cm	Violet	LV 0272	
	4	1	1 x 75 cm	Violet	LV 1273	
	5	2	1 x 75 cm	Violet	LV 1272	
	2	3/0	1 x 90 cm	Violet	LV 0373	
	3	2/0	1 x 90 cm	Violet	LV 0374	
	3.5	0	1 x 90 cm	Violet	LV 0375	
	4	1	1 x 90 cm	Violet	LV 0376	
HR 37 	3.5	0	1 x 75 cm	Violet	LV 0579	
	4	1	1 x 75 cm	Violet	LV 0379	
	3.5	0	1 x 90 cm	Violet	LV 0378	









CHIRASORB braided

UK04
70

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
 <p>HR 40</p>	4	1	1 x 45 cm	Violet	LV 0177	
	5	2	5 x 45 cm	Violet	LV 2250	
	3	2/0	1 x 75 cm	Violet	LV 0275	
	3.5	0	1 x 75 cm	Violet	LV 0276	
	4	1	1 x 75 cm	Violet	LV 0277	
 <p>HR 43</p>	3	2/0	1 x 90 cm	Violet	LV 1475	
	3.5	0	1 x 90 cm	Violet	LV 1476	
	4	1	1 x 90 cm	Violet	LV 1477	
	5	2	1 x 90 cm	Violet	LV 1478	
	3	2/0	1 x 120 cm	Violet	LV 1479	
 <p>HR 48</p>	4	1	1 x 75 cm	Violet	LV 1396	
	5	2	1 x 75 cm	Violet	LV 1400	
	3.5	0	1 x 75 cm	Violet	LV 1402	
 <p>HR 48 A</p>	3.5	0	1 x 90 cm	Violet	LV 1403	
	4	1	1 x 90 cm	Violet	LV 1401	
	5	2	1 x 90 cm	Violet	LV 1399	
 <p>HR 48 A</p>	4	1	4 x 45 cm	Violet	LV 0793	Take-off needle

Absorbable - Chirasorb braided






CHIRASORB braided

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 50 	5	2	5 x 45 cm	Violet	LV 1201	
	3.5	0	1 x 90 cm	Violet	LV 1392	
HR 60 	3 3.5	2/0 0	1 x 75 cm 1 x 75 cm	Violet Violet	LV 0278 LV 0279	
HR 65 	3 4 5	2/0 1 2	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	LV 1283 LV 1285 LV 1286	
	5	2	1 x 150 cm	Violet	LV 1287	
HR 75 	5	2	1 x 75 cm	Violet	LV 0496	
HR s: ½ circle, round body needle - slim						
HR s 18 	1.5	4/0	1 x 75 cm	Violet	LV 1256	Slim needle
HR s 27 Black 	3	2/0	1 x 75 cm	Violet	LV 2267	Slim, black needle



CHIRASORB braided

UK04
72





Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
HRX 22 	2	3/0	1 x 75 cm	Violet	LV 0420	
HRX 26 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	LV 0419 LV 0417	
HRX 27 	2 3 4	3/0 2/0 1	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	LV 0429 LV 0430 LV 0432	
HRX 30 	2	3/0	1 x 75 cm	Violet	LV 0426	
HRX 35 	3 4	2/0 1	1 x 75 cm 1 x 75 cm	Violet Violet	LV 0280 LV 0282	
	3.5	0	1 x 90 cm	Violet	LV 0381	

Absorbable - Chirasorb braided



Chirasorb braided

CHIRASORB braided










Needle	EP	USP	Length	Colour	Cat. No	Other information
HRX: ½ circle, round body needle - extra strong						
HRX 36 	3.5	0	1 x 75 cm	Violet	LV 1281	
	4	1	1 x 75 cm	Violet	LV 1282	
	3	2/0	1 x 90 cm	Violet	LV 1380	
HRX 40 	3.5	0	1 x 75 cm	Violet	LV 0285	
	4	1	1 x 75 cm	Violet	LV 0286	
	3.5	0	1 x 90 cm	Violet	LV 0385	
	4	1	1 x 90 cm	Violet	LV 0386	
	5	2	1 x 90 cm	Violet	LV 0387	
	3.5	0	1 x 150 cm	Violet	LV 4035	Loop
	4	1	1 x 150 cm	Violet	LV 4040	Loop
HRX 45 	4	1	1 x 90 cm	Violet	LV 0390	
	5	2	1 x 90 cm	Violet	LV 0391	
HRX 50 	3.5	0	1 x 75 cm	Violet	LV 0294	
	4	1	1 x 75 cm	Violet	LV 0295	
	5	2	1 x 75 cm	Violet	LV 0296	
	4	1	1 x 90 cm	Violet	LV 0395	
	5	2	1 x 90 cm	Violet	LV 0396	
	5	2	1 x 110 cm	Violet	LV 0393	











CHIRASORB braided

UK04
74

Absorbable - Chirasorb braided

Needle	EP	USP	Length	Colour	Cat. No	Other information
HRT: ½ circle, round body needle - tapercut						
HRT 17 	1.5	4/0	1 x 45 cm	Violet	LV 0566	
	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Violet Violet	LV 0561 LV 0562	
HRT 20 	1	5/0	1 x 45 cm	Violet	LV 0568	
HRT 22 	2	3/0	1 x 75 cm	Violet	LV 0572	
HRT 25 	2	3/0	1 x 75 cm	Violet	LV 0581	
HRT 31 	3	2/0	1 x 75 cm	Violet	LV 0606	
HRT 37 	3 3.5	2/0 0	1 x 75 cm 1 x 75 cm	Violet Violet	LV 0592 LV 0593	
	3 3.5 4	2/0 0 1	1 x 90 cm 1 x 90 cm 1 x 90 cm	Violet Violet Violet	LV 0597 LV 0598 LV 0599	
HRT 40 	3.5	0	1 x 75 cm	Violet	LV 0653	
HRT 45 	5	2	1 x 75 cm	Violet	LV 0665	
HRN: ½ circle, round body needle - blunt point						
HRN 50 	4 5	1 2	1 x 90 cm 1 x 90 cm	Violet Violet	LV 0686 LV 0675	

CHIRASORB braided







Needle	EP	USP	Length	Colour	Cat. No	Other information
HRF: ½ circle, flattened body needle, taper point						
HRF 17 	1 1.5	5/0 4/0	1 x 75 cm 1 x 75 cm	Violet Violet	LV 9242 LV 9256	
HRF 20 	1.5	4/0	1 x 75 cm	Violet	LV 9261	
HRF 26 	2 3	3/0 2/0	1 x 75 cm 1 x 75 cm	Violet Violet	LV 9266 LV 9267	
HRF 31 	3	2/0	1 x 75 cm	Violet	LV 9273	
DR: 3/8 circle, round body needle						
DR 15 	1.5	4/0	1 x 75 cm	Violet	LV 0724	
DR 21 	3	2/0	1 x 75 cm	Violet	LV 0744	
DR 30 	2	3/0	1 x 75 cm	Violet	LV 0746	
DRN: 3/8 circle, round body needle - blunt point						
DRN 60 	4	1	1 x 100 cm	Violet	LV 0699	













CHIRASORB braided

UK04
76

Absorbable - Chirasorb braided

Needle	EP	USP	Length	Colour	Cat. No	Other information
FR: 5/8 circle, round body needle						
FR 27	2	3/0	1 x 75 cm	Violet	LV 0609	
	3	2/0	1 x 75 cm	Violet	LV 0610	
	3.5	0	1 x 75 cm	Violet	LV 0611	
GR: Straight, round body needle						
GR 20	1.5	4/0	1 x 75 cm	Violet	LV 0402	
	3	2/0	1 x 90 cm	Violet	LV 0399	
KR: Ski, round body needle						
KR 25	3	2/0	1 x 20 cm	Violet	LV 2736	
						
KR 27	2	3/0	1 x 40 cm	Violet	LV 4630	Loop
	4	1	1 x 120 cm	Violet	LV 2390	
JRT: Hook shaped, round body needle - tapercut						
JRT 30	4	1	1 x 75 cm	Violet	LV 0694	
						
J-RX: Hook shaped, round body needle						
J-RX 30	3	2/0	1 x 75 cm	Violet	LV 1387	
						














CHIRASORB braided

Needle	EP	USP	Length	Colour	Cat. No	Other information
HS: ½ circle, reverse cutting needle						
 HS 18	2	3/0	1 x 75 cm	Violet	LV 0424	
	3	2/0	1 x 75 cm	Violet	LV 0425	
 HS 21	1.5	4/0	1 x 75 cm	Violet	LV 0300	
	2	3/0	1 x 75 cm	Violet	LV 0301	
 HS 25	2	3/0	1 x 75 cm	Violet	LV 0441	
	3	2/0	1 x 75 cm	Violet	LV 0442	
	3.5	0	1 x 75 cm	Violet	LV 0443	
 HS 26	3.5	0	1 x 75 cm	Violet	LV 1455	
 HS 30	3	2/0	1 x 75 cm	Violet	LV 0445	
	3.5	0	1 x 75 cm	Violet	LV 0446	
	4	1	1 x 75 cm	Violet	LV 0447	
 HS 36	3	2/0	1 x 75 cm	Violet	LV 0350	
	3.5	0	1 x 75 cm	Violet	LV 0351	
 HS 37	3.5	0	1 x 75 cm	Violet	LV 1451	
	4	1	1 x 75 cm	Violet	LV 1452	
	5	2	1 x 90 cm	Violet	LV 0453	
 HS 40	3	2/0	1 x 75 cm	Violet	LV 0455	
	3.5	0	1 x 75 cm	Violet	LV 0456	
	4	1	1 x 75 cm	Violet	LV 0457	
	5	2	1 x 75 cm	Violet	LV 0458	
	3.5	0	1 x 90 cm	Violet	LV 0460	
	4	1	1 x 90 cm	Violet	LV 0461	
 HS 48	3.5	0	1 x 75 cm	Violet	LV 0466	
	4	1	1 x 90 cm	Violet	LV 0469	
	5	2	1 x 90 cm	Violet	LV 0470	







CHIRASORB braided

UK04
78

Needle	EP	USP	Length	Colour	Cat. No	Other information
HSX: 1/2 circle, reverse cutting needle - extra strong						
HSX 26 	4 5	1 2	1 x 75 cm 1 x 75 cm	Violet Violet	LV 1456 LV 1454	
HSX 48 	3.5	0	1 x 90 cm	Violet	LV 0464	
DS: 3/8 circle, reverse cutting needle						
DS 11 	0.5	7/0	1 x 45 cm	Violet	LV 0193	
DS 12 	0.7 1 1.5	6/0 5/0 4/0	1 x 45 cm 1 x 45 cm 1 x 45 cm	Violet Violet Violet	LV 0201 LV 0202 LV 0203	
2 x DS 12 	1.5	4/0	1 x 45 cm	Violet	LV 1210	
DS 15 	0.7 1 1.5 2	6/0 5/0 4/0 3/0	1 x 45 cm 1 x 45 cm 1 x 45 cm 1 x 45 cm	Violet Violet Violet Violet	LV 0200 LV 0205 LV 0206 LV 0207	
DS 15 	1.5	4/0	1 x 75 cm	Violet	LV 0209	
DS 19 	1 1.5 2	5/0 4/0 3/0	1 x 45 cm 1 x 45 cm 1 x 45 cm	Violet Violet Violet	LV 0214 LV 0215 LV 0216	
DS 19 	1.5 2 3	4/0 3/0 2/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	LV 0211 LV 0212 LV 0213	
DS 24 	2	3/0	1 x 75 cm	Violet	LV 1225	
DS 25 	1.5 2 3	4/0 3/0 2/0	1 x 75 cm 1 x 75 cm 1 x 75 cm	Violet Violet Violet	LV 0222 LV 0223 LV 0224	
DS 25 	2 3	3/0 2/0	1 x 90 cm 1 x 90 cm	Violet Violet	LV 2223 LV 2224	
DS 26 	2	3/0	1 x 45 cm	Violet	LV 1221	

Absorbable - Chirasorb braided

CHIRASORB braided

Needle	EP	USP	Length	Colour	Cat. No	Other information
DS: 3/8 circle, reverse cutting needle						
 DS 30	3	2/0	1 x 75 cm	Violet	LV 0234	
	3.5	0	1 x 75 cm	Violet	LV 0235	
	3	2/0	1 x 90 cm	Violet	LV 1239	
	3.5	0	1 x 90 cm	Violet	LV 1235	
	4	1	1 x 90 cm	Violet	LV 1236	
DSC: 3/8 circle, conventional cutting needle						
 DSC 16	1.5	4/0	1 x 75 cm	Violet	LV 8208	
GS: Straight, reverse cutting needle						
 GS 13	1.5	4/0	1 x 45 cm	Violet	LV 1198	
Minimal invasive surgery						
	3	2/0	1 x 50 cm	Violet	MS 0143	Loop
	3.5	0	1 x 50 cm	Violet	MS 0144	Loop
 KR 27	4	1	1 x 120 cm	Violet	MS 1125	Loop



CHIRASORB braided natural ligatures

UK04
80

Pre-cut threads

EP	USP	Length	Cat. No
1	5/0	3 x 45 cm	LN 0012
1.5	4/0	3 x 45 cm	LN 0013
2	3/0	3 x 45 cm	LN 0014
3	2/0	3 x 45 cm	LN 0015
3.5	0	3 x 45 cm	LN 0016
2	3/0	10 x 45 cm	LN 0064
3	2/0	10 x 45 cm	LN 0065
1.5	4/0	2 x 75 cm	LN 0023
2	3/0	2 x 75 cm	LN 0024
3	2/0	2 x 75 cm	LN 0025



Chirapack - 20 pcs./box

EP	USP	Length	Cat. No
1.5	4/0	1 x 300 cm	LN 0003
2	3/0	1 x 300 cm	LN 0004
3	2/0	1 x 300 cm	LN 0005
3.5	0	1 x 300 cm	LN 0006













Absorbable - Chirasorb braided

CHIRASORB braided

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
HR 13 	1	5/0	1 x 45 cm	Natural	LN 0252	
	1.5	4/0	1 x 45 cm	Natural	LN 0253	
<hr/>						
HR 18 	1.5	4/0	1 x 75 cm	Natural	LN 0153	
	2	3/0	1 x 75 cm	Natural	LN 0256	LN 0257
<hr/>						
HR 22 	2	3/0	1 x 75 cm	Natural	LN 0262	
	3	2/0	1 x 75 cm	Natural	LN 0263	
<hr/>						
HR 27 	2	3/0	1 x 75 cm	Natural	LN 0266	
<hr/>						
HR 30 	3.5	0	1 x 75 cm	Natural	LN 0274	

CHIRASORB braided





Needle	EP	USP	Length	Colour	Cat. No	Other information
HRT: ½ circle, round body needle - tapercut						
HRT 17 	1.5	4/0	1 x 75 cm	Natural	LN 0561	
	2	3/0	1 x 75 cm	Natural	LN 0562	
HRT 22 	2	3/0	1 x 75 cm	Natural	LN 0572	
DS: 3/8 circle, reverse cutting needle						
DS 12 	0.7	6/0	1 x 45 cm	Natural	LN 1201	
	1	5/0	1 x 45 cm	Natural	LN 1202	
DS 15 	1.5	4/0	1 x 45 cm	Natural	LN 1206	
	2	3/0	1 x 45 cm	Natural	LN 1207	
DS 16 	1	5/0	1 x 45 cm	Natural	LN 1208	
DS 19 	1	5/0	1 x 45 cm	Natural	LN 1215	
	1.5	4/0	1 x 45 cm	Natural	LN 1216	
	2	3/0	1 x 45 cm	Natural	LN 1217	
DS 19 	1.5	4/0	1 x 75 cm	Natural	LN 1211	
	2	3/0	1 x 75 cm	Natural	LN 1212	
	3	2/0	1 x 75 cm	Natural	LN 1213	
DS 25 	2	3/0	1 x 45 cm	Natural	LN 1332	
	1.5	4/0	1 x 75 cm	Natural	LN 1231	
DS 25 	2	3/0	1 x 75 cm	Natural	LN 1232	
	3	2/0	1 x 75 cm	Natural	LN 1233	
DS 30 	3	2/0	1 x 75 cm	Natural	LN 1243	
	3.5	0	1 x 75 cm	Natural	LN 1244	



CHIRASORB braided

UK04
82

Absorbable - Chirasorb braided

Needle	EP	USP	Length	Colour	Cat. No	Other information
DSM: 3/8 circle, reverse cutting needle - micropoint						
DSM 12 	0.7	6/0	1 x 45 cm	Natural	LN 0201	Advanced cutting needle
	1	5/0	1 x 45 cm	Natural	LN 0202	
	1.5	4/0	1 x 45 cm	Natural	LN 0203	
DSM 15 	1	5/0	1 x 45 cm	Natural	LN 0206	Advanced cutting needle
	1.5	4/0	1 x 45 cm	Natural	LN 0207	
	2	3/0	1 x 45 cm	Natural	LN 0208	
DSM 19 	1.5	4/0	1 x 45 cm	Natural	LN 0216	Advanced cutting needle
	1.5	4/0	1 x 75 cm	Natural	LN 0212	
	2	3/0	1 x 75 cm	Natural	LN 0213	
DSM 25 	2	3/0	1 x 75 cm	Natural	LN 0223	Advanced cutting needle
	3	2/0	1 x 75 cm	Natural	LN 0224	