



Chiraplus™ Infusion set

Chiraplus Infusion sets include 15 µm filter, ventilation port and Luer lock adapter. All sets have 150 cm tubing, unless specified in the table. Chiraplus P is intended for pressure infusion and does not include flash ball. Chiraplus G is intended for gravity infusion and includes flash ball.

Product	Chamber Vol. (ml)	Packaging Type	Pcs./ Carton	Article No.
Chiraplus P	10.5	Blister	100	CHP336
Chiraplus P Dehp free	10.5	Blister	100	CHP536
Chiraplus P 12	12	Blister	100	CHP306
Chiraplus P 12 Dehp free	12	Blister	100	CHP506
Chiraplus P 190	10.5	Blister	100	CHP339
Chiraplus P 190 Safety	10.5	Blister	100	CHP338
Chiraplus P Dehp free Autofusion 180 *	15	Blister	300	CHP576
Chiraplus G	10.5	Blister	100	CHG331
Chiraplus G	10.5	Poly bag	500	CHG317
Chiraplus G Dehp free	10.5	Blister	100	CHG531
Chiraplus G 12	12	Blister	100	CHG301
Chiraplus G 12 Dehp free	12	Blister	100	CHG501
Chiraplus G	6	Poly bag	500	CHG318
Chiraplus G Dehp free	6	Poly bag	500	CHG518
Chiraplus G with needle	6	Poly bag	500	CHG319
Chiraplus G with metal spike	8.5	Poly bag	500	CHG316

Shelf life 5 years / *3 years

Chirahem™ Transfusion set

Chirahem Transfusion set is intended for blood transfusion administration. The set has a 200 µm blood filter.

Product	Chamber Vol. (ml)	Packaging Type	Pcs./ Carton	Article No.
Chirahem	17	Blister	100	CHT314
Chirahem Dehp free	17	Blister	100	CHT514
Chirahem	17	Poly bag	300	CHT320
Chirahem with needle	17	Poly bag	300	CHT321
Chirahem with metal spike	17	Poly bag	300	CHT322

Shelf life 5 years



According to EN, ISO standards. Quality management certified according to EN ISO 13485:2016. CE marked.