

Eye needles

Surgical Eye needles are made of high-quality stainless steel. Different types of needles are available to meet a full range of requirements. The needles are always packed in foil pouches with 12 pieces.

The surgical Eye needle range consists of Eye needles with a double spring eye and with a regular eye.

Eye needles
Double spring



Eye needles
Regular



Eye needles

Double spring



Needle	Needle no.	Wire diameter (mm)	Cat. No
PB HR: ½ circle, round body needle			
 <p>HR 61</p>	0	0.65	EY 8002
 <p>HR 52</p>	1	0.65	EY 8003
 <p>HR 42</p>	2	0.65	EY 8004
 <p>HR 36</p>	3	0.55	EY 8005
 <p>HR 30</p>	4	0.55	EY 8006



Eye needles

Double spring



UK04
204

Special Products - Eye needles

Needle	Needle no.	Wire diameter (mm)	Cat. No
PB HR: ½ circle, round body needle			
HR 24 	5	0.55	EY 8007
HR 20 	7	0.55	EY 8009
HR 13 	9	0.45	EY 8011
E HRX: ½ circle, round body needle - extra strong			
HRX 100 	0000	1.20	EY 8700
HRX 72 	000	1.10	EY 8701
HRX 63 	00	1.10	EY 8703
HRX 54 	0	1.00	EY 8705

Eye needles

Double spring



Needle	Needle no.	Wire diameter (mm)	Cat. No
E HRX: ½ circle, round body needle - extra strong			
HRX 46 	1	1.00	EY 8707
HRX 38 	2	0.85	EY 8709
HRX 30 	3	0.75	EY 8711
HRX 24 	4	0.65	EY 8713
HRX 22 	5	0.55	EY 8715
HRX 17.5 	6	0.55	EY 8717
PC GR: Straight, round body needle			
GR 70 	00	0.65	EY 8801
GR 52 	1	0.65	EY 8803
GR 36 	3	0.55	EY 8804
GR 24 	5	0.55	EY 8805

UK04
205



Special Products - Eye needles

Eye needles

Double spring



UK04
206

Needle	Needle no.	Wire diameter (mm)	Cat. No
GA HS: 1/2 circle, reverse cutting needle			

HS 72



0

1.80

EY 8623

HS 60



2

1.60

EY 8625

HS 54



3

1.40

EY 8626

HS 46



4

1.40

EY 8627

HS 38



5

1.30

EY 8628

HS 32



6

1.10

EY 8629

HS 23



8

1.00

EY 8631









Special Products - Eye needles

Eye needles

Double spring



Needle	Needle no.	Wire diameter (mm)	Cat. No
GE HS: ½ circle, reverse cutting needle			
<p>HS 80</p> 	1	1.30	EY 8601
<p>HS 68</p> 	3	1.10	EY 8603
<p>HS 58</p> 	5	1.10	EY 8605
<p>HS 50</p> 	8	1.00	EY 8608
<p>HS 42</p> 	11	0.85	EY 8611
<p>HS 34</p> 	14	0.75	EY 8614
<p>HS 28</p> 	17	0.65	EY 8617
<p>HS 24.5</p> 	19	0.65	EY 8619

UK04
207



Special Products - Eye needles

Eye needles

Double spring



UK04
208



Special Products - Eye needles

Needle	Needle no.	Wire diameter (mm)	Cat. No
GE HS: 1/2 circle, reverse cutting needle			
HS 20 	21	0.55	EY 8621
HS 17.5 	22	0.55	EY 8622
BE DS: 3/8 circle, reverse cutting needle			
DS 58 	5	1.10	EY 8405
DS 50 	8	1.00	EY 8408
DS 44 	10	0.85	EY 8410
DS 36 	13	0.75	EY 8413
DS 32 	15	0.75	EY 8415
DS 28 	17	0.65	EY 8417
DS 24.5 	19	0.65	EY 8419
DS 20 	21	0.55	EY 8421

Eye needles

Double spring







Needle	Needle no.	Wire diameter (mm)	Cat. No
BE DS: 3/8 circle, reverse cutting needle			
DS 17.5 	22	0.55	EY 8422
DS 12 	24	0.55	EY 8424

UK04
209

Eye needles

Regular



Needle	Needle no.	Wire diameter (mm)	Cat. No
PB HR: 1/2 circle, round body needle			
HR 24 	5	0.55	EY 9007
E HRX: 1/2 circle, round body needle - extra strong			
HRX 72 	000	1.10	EY 9701
HRX 38 	2	0.85	EY 9709
GA HS: 1/2 circle, reverse cutting needle			
HS 32 	6	1.10	EY 9629

Eye needles

Regular



UK04
210

Special Products - Eye needles

Needle	Needle no.	Wire diameter (mm)	Cat. No
GE HS: 1/2 circle, reverse cutting needle			
<p>HS 42</p>	11	0.85	EY 9611
<p>HS 34</p>	14	0.75	EY 9614
<p>HS 24.5</p>	19	0.65	EY 9619
<p>HS 20</p>	21	0.55	EY 9621
BE DS: 3/8 circle, reverse cutting needle			
<p>DS 58</p>	5	1.10	EY 9405
<p>DS 44</p>	10	0.85	EY 9410
<p>DS 32</p>	15	0.75	EY 9415
<p>DS 28</p>	17	0.65	EY 9417
<p>DS 24.5</p>	19	0.65	EY 9419
<p>DS 20</p>	21	0.55	EY 9421