



Vitrex® Alcohol Swabs

Vitrex Alcohol swabs are an easy and convenient way to cleanse the skin and help fight bacteria at the site of application. Use as an antiseptic skin cleanser before injections and blood sampling.

Each Alcohol Swab is individually sealed and disposable for convenience and safety.

Technical Specifications

Description	Size	Content	Pcs./Box	Article No.
Alcohol Swabs	30 x 60 mm	70% Isopropyl Alcohol	100	520213
Alcohol Swabs	30 x 60 mm	82% Ethanol	100	521213
Alcohol Swabs	30 x 60 mm	82% Ethanol + 0.5% Chlorhexidine	100	522213

Shelf life 5 years



Transit packaging

Pallet

Packaging:
 Primary packaging: 100 pcs.
 Transit packaging: 10,000 pcs.
 Pallet: 400,000 pcs.

Vitrex® Safety Containers

Vitrex Safety Containers are designed for hospitals and healthcare institutions for clinical disposal of lancets, needles, scalpels, syringes and other types of sharps.

Each type of lid is specially designed to help the user dispose needles, syringes, etc. Lid with reversible closure can be safely fixed.



Technical Specifications

Description	Size	Content	Article No.
Safety container	4.7 cm x 12.9 cm	0.2 l.	193800
Safety container	9 cm x 10.8 cm	0.5 l.	194000
Safety container	13.7 cm x 12.5 cm	1.5 l.	194200
Safety container	13.7 cm x 18 cm	2 l.	194300
Safety container	21.1 cm x 19.5 cm	5 l.	194700